## EN User manual for silicone foam dressing

### Description

Silicone foam dressing is an absorbent, atraumatic dressing made from polyurethane foam. The outer surface of the foam is bonded to a vapour-permeable polyurethane membrane, which acts as a barrier to liquid and microorganisms. The wound contact surface of Silicone foam dressing is coated with a layer of soft silicone that does not stick to the surface of a wound or cause trauma to delicate new tissue upon removal. This soft silicone layer is also slightly tacky, which facilitates application and retention of the dressing to intact skin, but does not cause epidermal stripping or pain on removal.

The centre of the dressing is polyurethane foam that can lock exudate into the product, preventing the wound bed from drying out and the surrounding skin from becoming macerated.

### Indications

The product is suitable for many types of exuding wounds, including leg and pressure ulcers, superficial and partial thickness burns, donor sites, postoperative wounds and skin abrasions. The dressing absorbs exudate and maintains a moist wound-healing environment whilst minimising the risk of maceration.

### Contra-indications

Discontinue use if the patient is allergic to any components of the dressing.

### Application method

1. Make sure the wound area is clean and dry before using the dressing.
2. A dressing should be selected that overlaps the wound margin by at least two centimetres.
3. Remove the release film before use and adhesive the dressing on the wound area.

### Frequency of dressing changes

The interval between changes will normally be determined by the degree of exudate produced. The dressing may be left up to 7 days.

### Warnings


Sloughy wounds dressed with silicone foam dressing may initially appear to increase in size due to autolytic debridement promoted by the moist conditions produced beneath the dressing. This is normal and to be expected. Not for use with oxidizing solutions, such as Hypochlorite solution or Hydrogen peroxide solution (H<sub>2</sub>O<sub>2</sub>).

### Storage

Silicone foam dressings are stored in a well-ventilated clean room with relative humidity less than 80%, and no corrosive gases.

### Other

 Silicone foam dressings have been sterilized using ethylene oxide.

-  Do not reuse or resterilize product.
-  Do not use if package has been compromised.
-  Protect against direct sunlight.
-  Keep dry.
-  Store at a normal room temperature 5 - 35 °C.
-  Instructions for use.

DECIFERA® is a registered trademark of Van Heek Medical. All rights reserved.

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