BBrain

**BASIC 10 INCH
CALENDAR CLOCK**

Manual

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Thank you for purchasing the BBrain clock. We hope that the clock contributes to facilitating the user's life.

1. About this manual


- Read all safety information and this manual carefully and make sure that you are handling the device safely and correctly before you take the BBrain into first use.
- The images and screenshots that are used in this manual may differ from the actual product.
- Available functions and extra services depend on the hardware and the software .
- BBrain is not responsible for performance or compatibility issues caused by modifications in the register made by the user.
- Audio sources, backgrounds and images within this device are provided through licenses with limited user rights between BBrain and the respective owners. Taking over and using these materials for commercial or other purposes is a violation of the copyright law. BBrain is not responsible for the violation of the copyright caused by the user.
- Store this manual carefully for future advise.

2. Description of the product

Check if all necessary parts are provided in the package before you continue with the installation of your BBrain clock. If this is not the case, please send an e-mail to: info@bbrain.eu

2.1 Content

- 1 BBrain calendar clock
- 1 BBrain stand
- 1 BBrain charger with charger cable (2 meter)
- 1 Manual

3. Startup

Connect the charger cable. Press the On/Off button at the back of the device. A startup screen will be showed for approximately 5 seconds, as showed in Figure 1.


Figure 1

After the startup screen, the default time appears as showed in Figure 2.


Figure 2

4. Installation of the calendar clock

Press MENU to enter the settings. The settings will be displayed in the screen. Start the navigation in the menu by pressing **⬆** and **⬇** and then choose the settings you want to adjust. When everything is set, click the MENU to return to the clock.

LANGUAGE	<ENGLISH>
-----------------	------------------------

Press **⬅** or **➡** to select the language: Dutch, German, French, Italian, Spanish, Polish, Danish, Norwegian, Finnish, Swedish and English. Press **⬇** to continue in the menu.

TIME	< 00:00 >
-------------	------------------------

Press **OK** to set the time. To mark the hours or minutes with a white underline press **⬅** or **➡**. 13:00, the first two digits left of the : are the hours. Use **⬆** or **⬇** set the hours. To set the minutes, press **⬅** or **➡** to select the minutes. 15:00 the last two digits of the : are the hours. Press the **⬆** or **⬇** to set the minutes. Press the **OK** button to confirm the time.

DATE	< 00-00-2019 >
-------------	-----------------------------

To set the date press **OK**. By pressing the buttons **⬅** or **➡** the day, month or year will be marked with a white underline. You can adjust each part by pressing **⬆** or **⬇**. Example: 01-01-2019 , use **⬆** or **⬇** to choose the right day. Then press the **⬅** or **➡** buttons to choose the month. The month will be marked 01-01-2019. Repeat the steps to change the year. Press **OK** to confirm the date.

TIME MODE	< 24 Hours >
------------------	---------------------------

Press the **OK** button to select the 12-hour (AM / PM) or 24-hour clock. It is important to choose the 24-hour mode before the right time is set, so that the clock 'knows' whether it is 3 o'clock in the afternoon or 3 o'clock at night. Is 3 hours PM 3 pm in the 24-hour mode. Press **⬇** to continue in the menu.

DATE MODE	< Month-Day-Year >
------------------	---------------------------------

Press **OK** to choose the following 3 views: day-month-year, month-day-year or year-month-day . Press **⏴** to continue in the menu.

THE MORNING STARTS AT	< 06 >
------------------------------	---------------------

To indicate the correct time for the start of the morning, press the **⏴** or **⏵**. Press **⏴** to continue in the menu.

THE AFTERNOON STARTS AT	< 12 >
--------------------------------	---------------------

To indicate the correct time for the start of the afternoon, press the **⏴** or **⏵**. Press **⏴** to continue in the menu.

THE EVENING STARTS AT	< 18 >
------------------------------	---------------------

To indicate the correct time for the start of the evening, press the **⏴** or **⏵**. Press **⏴** to continue in the menu.

THE EVENING NIGHT AT	< 00 >
-----------------------------	---------------------

To indicate the correct time for the start of the night, press the **⏴** or **⏵**. Press **⏴** to continue in the menu.

SELECT HOUR FOR TONE	< 1 2 3 4 5 6 7 8 9 10 11 12 13 14 15 16 17 18 19 20 21 22 23 00 >
-----------------------------	-------------------------------------------------------------------------------------

For every hour, there can be set a tone as an additional reminder. Press **⏴** or **⏵** if the hour is selected, the hour is marked with a white underline. Confirm the hour by pressing **OK** and continue in the menu. The programmed hours are displayed in blue. The memories can be deleted in the same way.

BRIGHTNESS DAY	< 07 >
-----------------------	---------------------

Set the daytime brightness. Press **⏴** or **⏵** to select the correct brightness. Press **⏴** to continue in the menu.

BRIGHTNESS NIGHT	< 01 >
-------------------------	---------------------

Set the nighttime brightness. Press **⏴** or **⏵** to select the correct brightness. Press **⏴** to continue in the menu.

VOLUME	< 04 >
---------------	---------------------

Press **⏪** or **⏩** to adjust the volume. Press **⏴** to continue in the menu.

COLOR MODE	< Black / Red Bar >
-------------------	----------------------------------

Press **⏪** or **⏩** to choose the following 3 views: Black / White Text (standard mode), Black / Red Bar (suitable for people with early dementia), or Black / Yellow Text (suitable for the visually impaired).

DISPLAY INTERVAL MODE	< OFF >
------------------------------	----------------------

Select the interval screen mode, choose 30 seconds, 60 seconds, 90 seconds or screen mode off. Press **⏴** to continue in the menu.

ALARM SETTING	< OFF >
----------------------	----------------------

Press **OK** to set various alarms. Press the again **OK** to set the time. Press **⏪** or **⏩** to mark the hours or minutes with a white underline. Use **⏴** or **⏵** to set the correct hours. Repeat to set the minutes. First press **⏪** or **⏩** to select the minutes. Then press **⏴** or **⏵** to set the minutes. Press **OK** to confirm the time. Choose the days to repeat the alarm. The scheduled days are shown in blue. Press **⏴** to select the duration: 30 seconds, 45 seconds, or 60 seconds. The steps above can be repeated for the other alarms.

- GO TO BED ALARM: See the description above
- TAKE A NAP ALARM: See the description above
- WAKE UP ALARM: See the description above
- BIRTHDAY USER ALARM - CHRISTMAS DAY ALARM - NEW YEAR'S DAY ALARM ON/OFF , day, month, hour and minutes can be set (only available at the SW4 and SW5 calendar clock)

MEDICINE SETTINGS	>
--------------------------	-------------

There are several medicine reminders that can be set. These will be explained further in this manual. Every reminder can be set to repeat 0, 1, 2 or 3 times. Also the duration can be set (30-60-90 seconds). When the arrow is selected and **OK** is clicked, the Medicine menu will be showed and can be set. For further explanation read the section 'Medicine Settings'. (only available at the SW4 and SW5 calendar clock)


RESET SETTINGS	>
-----------------------	-------------

When selected, press **OK** to reset all settings. Press **⏪** or **⏩** to select yes or no and confirm by pressing **OK**. If you are ready to adjust your settings, press MENU to exit the menu. The clock automatically closes the menu after 10 seconds.

SET DISPLAY	>
--------------------	---

To set the display, press **OK** and hold for 3 seconds.
The display shows the different views.


Only available at the analog SW5 clock.

EXAMPLES OF THE DIFFERENT ALARMS AND MEDICINE REMINDERS

When set the alarm(s) and medicine reminder(s), these will appear when the time is reached. These will be showed with a picture where the time is displayed with a message in this and a sound for the duration as set. If also a repeat is set, this will be repeated after five minutes.

Some examples of these messages will be showed below.

AM / PM are showed when the time mode is set to standard time, otherwise it will be military time. See Items in menu.

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- Geluidsbronnen, achtergronden en afbeeldingen in dit apparaat zijn onder licentie verstrekt met beperkte gebruiksrechten tussen BBrain en de betreffende eigenaren. Het overnemen en gebruiken van deze materialen voor commerciële of andere doeleinden maakt inbreuk op de copyright-wetgeving. BBrain is niet verantwoordelijk voor inbreuk op het copyright door de gebruiker.