



















MDR Symbology Reference Guide

Sutures







	Undyed, Monofilament, Non-Absorbable, Uncoated	ETHILON™, PROLENE™, Surgical Stainless Steel		Dyed, Braided, Non-Absorbable, Uncoated	MERSILENE™
	Undyed, Monofilament, Absorbable, Uncoated	MONOCRYL™, MONOCRYL™ Plus, PDS™ II, PDS™ Plus		Dyed, Braided, Absorbable, Uncoated	PDS™ Cord, VICRYL™
	Undyed, Braided, Non-Absorbable, Uncoated	MERSILENE™ Suture & Tape, Surgical Stainless Steel		Dyed, Braided, Non-Absorbable, Coated	ETHIBOND EXCEL™, PERMA-HAND™
	Undyed, Braided, Absorbable, Uncoated	VICRYL™		Dyed, Braided, Absorbable, Coated	Coated VICRYL™, Coated VICRYL™ Plus, VICRYL™ Rapide
	Undyed, Braided, Non-Absorbable, Coated	ETHIBOND EXCEL™, PERMA-HAND™		Dyed, Twisted, Non-Absorbable, Uncoated	Virgin Silk
	Undyed, Braided, Absorbable, Coated	Coated VICRYL™, Coated VICRYL™ Plus, VICRYL™ Rapide		Undyed, Barbed, Monofilament, Absorbable, Uncoated, Unidirectional	STRATAFIX™ Spiral MONOCRYL™ Plus
	Undyed, Twisted, Non-Absorbable, Coated	Temporary Cardiac Pacing Wire		Undyed, Barbed, Monofilament, Absorbable, Uncoated, Bidirectional	STRATAFIX™ Spiral MONOCRYL™ Plus
	Dyed, Monofilament, Non-Absorbable, Uncoated	ETHILON™, MERSILENE™, PROLENE™, PRONOVA™		Dyed, Barbed, Monofilament, Absorbable, Uncoated, Unidirectional	STRATAFIX™ Spiral MONOCRYL™ Plus, STRATAFIX™ Spiral PDS™ Plus, STRATAFIX™ Symmetric PDS™ Plus
	Dyed, Monofilament, Absorbable, Uncoated	MONOCRYL™, MONOCRYL™ Plus, PDS™ II, PDS™ Plus, VICRYL™		Dyed, Barbed, Monofilament, Absorbable, Uncoated, Bidirectional	STRATAFIX™ Spiral MONOCRYL™ Plus, STRATAFIX™ Spiral PDS™ Plus











Antibacterial

MONOCRYL™ Plus, VICRYL™ Plus, PDS™ Plus, STRATAFIX™ Spiral MONOCRYL™ Plus, STRATAFIX™ Spiral PDS™ Plus, STRATAFIX™ Symmetric PDS™ Plus

Mesh Devices

	Absorbable Mesh Device	VICRYL™ Knitted Mesh		Partially Absorbable Mesh Device	ARTISYN™ Y-Shaped Mesh, ULTRAPRO™ Mesh, ULTRAPRO ADVANCED™ Mesh, ULTRAPRO™ Hernia System, ULTRAPRO™ Plug
	Absorbable Adhesion Barrier	GYNECARE INTERCEED™		Permanent Mesh / Mesh Device with Absorbable Layer	PROCEED™ Surgical Mesh, PROCEED™ Ventral Patch
	Partially Absorbable Mesh / Mesh Device with Absorbable Layer	PHYSIOMESH™ Open		Permanent Mesh Device	GYNECARE GYNEMESH™ PS, PROLENE™ Mesh and Soft Mesh, PROLENE™ Hernia System, GYNECARE TVT™


BLAKE™ Drains and J-VAC™ Reservoirs

							
Suction Reservoir	Bulb Suction Reservoir	Drain Adaptor	Silicone Flat Drain	Silicone Round Drain	Cardio Connector 1:1	Cardio Connector 2:1	Cardio Connector 3:1













DERMABOND™ Topical Skin Adhesive

	
Contents: Vial	Crush Vial After Opening the Package

Adjunctive Hemostats

	Absorbable Hemostat	SURGICEL™ Powder, SURGICEL SNoW™, SURGICEL NU-KNIT™, SURGICEL FIBRILLAR™, SURGICEL™ Original
--	---------------------	--

ISO Standard Symbols

					
Packaging Unit	Date of Manufacture	MR Safe	Medical Device	Single sterile barrier system with protective packaging inside	Unique Device Identifier
					
Information website for patients	Patient name or patient ID	Date of implantation	Name/address of the implanting healthcare	Consult IFU or electronic IFU	Contains hazardous substance

Please refer always to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

Johnson & Johnson Medical GmbH
Robert-Koch-Str. 1, 22851
Norderstedt Germany
www.jnjmedtech.com

Not all products may be available in all Markets.
© Johnson & Johnson Medical GmbH, 2024,
EM_ETH_WOUN_392307, EMEA